

Cricoid Stick Trainer

Simulating realistic needle and surgical cricothyrotomy

Needle and surgical cricothyrotomy skills can be practiced on this model with interchangeable rigid and soft tracheas, that also articulates to a range of adult manikins.

261-01001 Cricoid Stick Trainer

Cricoid Stick Trainer Includes:

Head, Rigid Trachea with Simulated Lung, Soft Trachea with Simulated Lung, Replaceable Neck Skin, Base, and Directions for Use.



Consumables

261-01150
Replaceable
Trachea – Rigid



261-01250
Replaceable
Trachea – Soft



261-01350
Replacement
Neck Skin



Choking Charlie

Designed specifically for training of the Heimlich Abdominal Thrust Maneuver

The Choking Charlie manikin is an adult torso designed specifically for training students in the performance of the Heimlich Abdominal Thrust Maneuver. Realistic anatomy and response using simulated boluses provide instructors with an excellent tool for instructing and practicing this lifesaving technique.

102-00001 Choking Charlie

Choking Charlie Manikin Includes:

Torso, 4 Simulated Boluses, talcum powder, Tank Top, Carry Case, and Directions for Use.



Consumables

102-00150
Simulated
Boluses (4)



Pneumothorax Trainer



The Pneumothorax Trainer, with bilateral tension pneumothoraces and decompression portals in each axilla and subclavian region, is specifically designed for training professionals in the practice of chest decompression.

260-05001 Pneumothorax Trainer

Pneumothorax Trainer Includes:

Torso Manikin, Bladder Replacement Kit, Airway Lubricant, Tank Top, Carry Case, and Directions for Use.

Consumables

250-21050
Airway
Lubricant
(45ml)



200-00150
Bladder Replacement Kit



Mr. Hurt Head Trauma Trainer

Realistic cranial trauma management

Cast from an actual victim of a highway accident, Mr. Hurt Head Trauma Trainer presents students with a variety of trauma conditions.

This trainer comes mounted to a base, but can be easily transferred to adult manikins for use in full-body trauma scenarios.



261-02001 Mr. Hurt Head Trauma Trainer