

# Resusci Anne Advanced SkillTrainer

Advanced skills and quality CPR training



SimPad is included



See pages 32-33 for more information regarding the SimPad® System and Patient Monitor. Sold separately.

The Resusci Anne Advanced SkillTrainer offers high-quality, yet relatively low-cost, advanced skills training. The combination of the remote control and PC SkillReporting software allows comprehensive instructor-to-student debriefing and documentation of key events. Additionally, first aid and extrication practice can be incorporated into the learning by adding the optional first aid and rescue limbs.

151-000xx Resusci Anne Advanced SkillTrainer

### Resusci Anne Advanced SkillTrainer Includes:

Full-Body Manikin with Track Suit, Carry Case, IV Arm, Software CD, USB Interface Cable, 1 Bottle Simulated Blood, Airway Lubricant, SimPad and Directions for Use.

xx = Language versions

## Consumables

250-21050  
Airway  
Lubricant  
(45 ml)



300-00750  
Simulated  
Blood



325-00350  
Replacement Skin and Vein  
System – Female

